

CARTIVA[®]

Synthetic Cartilage Implant

The Only **PMA Approved** Product for
the Treatment of 1st MTP Osteoarthritis




THE DIFFERENCE IS MOVING.™

CARTIVA[®]

THE DIFFERENCE IS MOVING.™

The Only **PMA Approved** Product for
the Treatment of 1st MTP Osteoarthritis




THE FIRST & ONLY PMA ALTERNATIVE TO FUSION + IMPROVED RANGE OF MOTION



SUBSTANTIAL REDUCTION IN PAIN

A substantial and clinically meaningful reduction in pain using the Visual Analog Scale (VAS) was experienced by subjects in the Cartiva implant group at every follow-up visit through 2 years. Cartiva implant subjects demonstrated a 93% reduction from a median score of 68 at baseline to 5 at 2 years.



SUBSTANTIAL FUNCTIONAL IMPROVEMENT

Functional activities were evaluated using the validated Foot and Ankle Mobility Measure (FAAM). Substantial improvement was observed for the Cartiva implant subjects throughout the 2-year follow-up period with a 168% median improvement observed in the sporting activities scale.



IMPROVED RANGE OF MOTION



There was a substantial and clinically important improvement in median active dorsiflexion motion in the Cartiva implant group, restoring motion to levels which are documented in the literature to be needed for normal walking gait while experiencing substantial reduction in pain.

Level I Clinical Evidence¹ of safety and effectiveness for treatment of 1st MTP Osteoarthritis, in the largest randomized study ever conducted for this condition.



EXCELLENT SURVIVORSHIP: **96.2%** implants retained

HIGH PATIENT SATISFACTION: **96%** of patients would undergo procedure again



FASTER THAN FUSION

FAST & SIMPLE SURGICAL PROCEDURE


Cartiva® surgeries are **40%** (23 minutes) faster than fusion surgeries.



Illustration of the Cartiva device implanted into metatarsal head



Damaged cartilage replaced with new Cartiva implant bearing surface



In most operating rooms in the United States, the value of a minute can be as high as \$100.²

PATIENT BENEFITS

QUICKER RECOVERY

Cartiva® SCI patients return to pre-operative activities faster than fusion patients.



- No cast, full weight bearing immediately as tolerated, able to drive
- Range of motion exercises encouraged immediately

HYDROGEL THAT WORKS LIKE NATURAL CARTILAGE


Mechanical and physical properties similar to native cartilage.

PROPERTY	ARTICULAR CARTILAGE ^{1,3}	CARTIVA
Water Content	60-80%	60%
Compressive Modulus	0.3 – 0.8 MPa	2.5-3.2 MPa
Coefficient of Friction	<0.01 – 0.05	0.04 – 0.07

FEATURES	BENEFITS
Synthetic	No risk of viral or bacterial transmission associated with human or animal derived materials
Biocompatible	Composed of saline and an organic polymer
Durable	Mechanical and physical properties similar to native cartilage capable of withstanding repetitive loading typical of MTP joint
Slippery	Low coefficient of friction aids joint articulation and mobility

Patient Satisfaction

% of Patients that **would** have the procedure again.



PROVEN RESULTS



CLINICAL STUDIES

Patients experience substantial reduction in pain, function improvement, and increased range of motion.

2 YEARS



N = 130

Substantial Pain Reduction




REDUCTION

Substantial Functional Improvement



IMPROVEMENT

Improved Range of Motion



IMPROVEMENT

EXTENSIVELY TESTED

BIOCOMPATIBILITY OF CARTIVA DEVICE					
Test	Method/Model		Result		
Cytotoxicity	L929 MEM Elution		Non-cytotoxic		
Cytotoxicity	Direct Contact		Non-cytotoxic		
Sensitization	Kligman Maximization		Non-sensitizer		
Irritation/Intracutaneous	IC Injection		Negligible irritant		
Acute Systemic Toxicity	Systemic Injection		Negative		
Subchronic Toxicity	Femoral Condyle Implantation		Non-toxic		
Chronic Toxicity	Femoral Condyle Implantation		Non-toxic		
Genotoxicity	Ames Reverse Mutation		Non-mutagenic		
Genotoxicity	Chromosomal Aberration Assay		Non-clastogenic		
Genotoxicity	Rodent Bone Marrow Micronucleus		Non-clastogenic		
Implantation	Bone Implantation In Femoral Condyle		Negative/no reaction		
Pyrogenicity	Rabbit Pyrogen Test		Non-pyrogenic		
BIOCOMPATIBILITY OF CARTIVA INSTRUMENTATION					
Cytotoxicity	L929 MEM Elution		Non-cytotoxic		
Sensitization	Kligman Maximization		Non-sensitizer		
Irritation/Intracutaneous	IC Injection		Negligible irritant		
ANIMAL SAFETY STUDIES					
Animal Study 1 Year Goat	Cartiva device implanted in load bearing region of medial femoral condyle in stifle of 8 mature goats; control defects in 4 goats At one year, knees evaluated via - High field strength MR imaging system for morphology and quantitative T2 and T1-rho parameters; - Histological processing - Biomechanical testing		<ul style="list-style-type: none"> - No evidence of local or systemic toxicity - No inflammatory reaction around implant or osteolytic bone loss - Non-significant change to opposing tibial-surface - No difference in presence of subarticular cysts with control - No device fragmentation or dislodgement - No particulate migration 		
Particulate Implant Study 6 month rabbit	<ul style="list-style-type: none"> - 5 million cycle wear debris quantified and characterized - Particulate replicated and injected via bolus in a quantity 9x - Test injections and control (saline) administered to 16 animals. At 3 and 6 months, histology and pathology per ISO standards		<ul style="list-style-type: none"> - No complications on injection - No test-article related adverse changes - No significant findings on clinical observation, gross pathology, histomorphometry, or histopathology of localized tissue - Systemic issues showed no microscopic changes related to the treatment - No wear debris or foreign body giant cells with injected material 		
FUNCTIONAL TESTING					
Fatigue Testing	Cycles Test Surface Axial Load	5 million Stainless Steel 4 MPa	<ul style="list-style-type: none"> - Mechanical durability demonstrated after 5M continuous cycles at peak load of 4 MPa - Significant mass and height recovery upon unloading - The Cartiva device demonstrated adequate strength to survive the repetitive, compressive loads that occur clinically in the 1st MTP. 		
Wear Testing	Cycles Test Surface Simulated Axial Load	5 million Cartilage 4 MPa	<ul style="list-style-type: none"> - Resistance to wear demonstrated after 5M continuous cycles at simulated peak load of 4 MPa - 0.18% average mass loss (1.64mg) - Worse case wear debris over 5 years of 2.88 mg or 0.31% - Volumetric wear rate of 1.50mm³/yr that is considerably lower than UHMWPE (80mm³/year)⁴ 		
MATERIALS PROPERTIES					
Unconfined Compression	Loading of unconfined devices to achieve 10%, 20%, 30% and 40% strain to measure deformation resistance of the matrix and determine compatibility of the device with surrounding native tissue			CARTIVA	Articular Cartilage
			Compressive Modulus	3.05 ± 0.12 MPa	.31–.80 ⁵ MPa
			Equilibrium Elastic Compressive Modulus	2.68–3.34 MPa	0.54 ⁶ MPa
Confined Compression	Devices confined in compression fixture with 5%, 10%, 15%, 20% and 25% strain applied to assess matrix stiffness at equilibrium (ie when load-induced fluid flow has ceased).		Higher polymer content and presence of physical cross links in Cartiva results in a mean aggregate modulus of 6.7 ± 1.0 MPa where cartilage values range between 0.6 and 1.2 MPa.		
Shear	Devices seated between test blocks that are moved apart perpendicularly until failure or 5 mm displacement; thereby, providing a baseline understanding of the simple shear properties of the material.			CARTIVA	Articular Cartilage
			Shear Moduli	0.16–0.36 MPa	0.45 ⁷ MPa (0.22–0.68 MPa)
			Fatigued devices exhibited no change in shear properties and resistance to mechanically induced degradation properties. All devices exhibited full 100% lateral shear strain without tearing or showing shear fracture.		
Creep	4 MPa loading in confined compression fixture to elucidate structural changes since equilibrium swelling properties are sensitive to the nature and stability of the hydrogel crosslinks		<ul style="list-style-type: none"> - Biphasic creep - 4-5% mass loss 		
S-N Analysis	Devices loaded in a confined fixture to 8, 12, 18, and 24 MPa out to 5,000,000 cycles		<ul style="list-style-type: none"> - No catastrophic failure - Continuous 5M compression cycles - Extreme loads of 24 MPa (6 x peak load) - Even under significant stresses, no failures 		

ORDERING INFORMATION

For Customer Service Call: 877-336-4616

IMPLANTS

CAR-10-US

10 mm Cartiva MTP Implant

**CAR-8-US**

8 mm Cartiva MTP Implant



DRILL BITS

MTD-10

10 mm Counterbore Drill Bit

**MTD-8**

8 mm Counterbore Drill Bit



GUIDE PINS

PNN-02

2 mm Guide Pin, Non-Threaded
(6 per pack)



INTRODUCERS

INT-10

10 mm Introducer

**INT-8**

8 mm Introducer



PLACERS

PLC-10

10 mm Placer

**PLC-8**


8 mm Placer



DELIVERY TRAY

TRA-05-US

Delivery Tray



Brief Summary of Important Product Information

INDICATIONS

The Cartiva® Synthetic Cartilage Implant is intended for use in the treatment of patients with painful degenerative or post-traumatic arthritis (hallux limitus or hallux rigidus) in the first metatarsophalangeal joint with or without the presence of mild hallux valgus, defined as a hallux valgus angle less than or equal to 20° (greater than 20° was an exclusion criteria in the clinical study).

CONTRAINDICATIONS

The Cartiva SCI should not be implanted in subjects with the following conditions:

- Active infection of the foot
- Known allergy to polyvinyl alcohol
- Inadequate bone stock due to significant bone loss, avascular necrosis, and/or large osteochondral cyst (> 1 cm) of the metatarsophalangeal joint
- Lesions of the first metatarsal head greater than 10 mm in size
- Diagnosis of gout with tophi
- Physical conditions that would tend to eliminate adequate implant support (e.g., insufficient quality or quantity of bone resulting from cancer, congenital dislocation, or osteoporosis), systemic and metabolic disorders leading to progressive deterioration of bone (e.g., cortisone therapies or immunosuppressive therapies), and/or tumors of the supporting bone structures

PRECAUTIONS

The safety and effectiveness of this device has not been established in subjects with the following conditions:

- Pediatric patients (< 22 years of age)
- Subjects with osteonecrosis of the first metatarsophalangeal joint
- Osteoarthritis involving the first metatarsophalangeal joint with grade 0 or 1 hallux rigidus per the Coughlin Scale^a

CITATIONS:

1. Data on file at Cartiva, Inc.
2. Daniels TR, Younger ASE, Penner MJ, Wing MJ, Miniaci-Coxhead SL, Pinsky E, Glazebrook M. Mid-term Outcomes of Polyvinyl Alcohol Hydrogel Hemiarthroplasty of the First Metatarsophalangeal Joint in Advanced Hallux Rigidus. Foot Ankle Int. First published on November 30, 2016
3. Garner, Patrick. Complexities in the Operating Room. Industrial and Systems Engineering Research Conference. Proceedings (2012):1-8. Web address accessed December 6, 2016: <https://www.iienet2.org/uploadedfiles/SEMS/Students/GarnerPITT2012.pdf>
4. Baker MI, Walsh SP, Zvi Sc, Boyan BD, A Review of polyvinyl alcohol and its uses in cartilage and orthopedic application, J Biomed Mater Res B Appl. Biomater. 2012 Jul; 100(5):1451-7
5. Jacobs CA, Christensen CP, Greenwald AS, McKellop H, Clinical performance of highly cross-linked polyethylenes in total hip arthroplasty. J Bone Joint Surg Am, 2007;89(12):2779-2786
6. Korhonen RK, Laasanen MS, Toyras J, Rieppo J, Hirvonen J, Helminen JF, Jurvelin JS, Comparison of the Equilibrium Response of Articular Cartilage in Unconfined Compression, Confined Compression and Indentation, J Biomech. 2002 Jul;35(7):903-909
7. Jurvelin JS, Buschmann MD, Hunziker EB, Optical and Mechanical Determination of Poisson's Ratio, J Biomechanics. 1997;30(3):235-241
8. Athanasiou KA, Liu GT, Lavery LA, Lanctot DR, Schenck RC Jr, Biomechanical Topography of Human Articular Cartilage in the First Metatarsophalangeal Joint, Clin Orthop Relat Res. 1998 Mar;(348):269-281
9. Coughlin MJ, Shurnas PS. Hallux rigidus. Grading and long-term results of operative treatment. American Journal of Bone Joint Surgery. 85-A(11):2072-88. November 2003

CARTIVA®

Cartiva, Inc.
6120 Windward Parkway, Suite 220
Alpharetta, GA 30005

(877)336-4616

www.cartiva.net
info@cartiva.net

© 2017 Cartiva, Inc. All rights reserved. Patent:: <http://cartiva.net/Home/Patents>

B20-0368 Rev. E

The safety and effectiveness of the Cartiva SCI device for treatment in the presence of hallux varus to any degree or hallux valgus >20° is unknown.

The safety and effectiveness of using more than one Cartiva SCI device per joint is unknown.

The safety and effectiveness of the Cartiva SCI device at anatomic locations other than the first metatarsophalangeal joint is unknown.

The Cartiva SCI device should only be used by experienced surgeons who have undergone training in the use of this device. A lack of adequate experience and/or training may lead to a higher incidence of adverse events.

Examine all instruments prior to surgery for wear or damage. Replace any worn or damaged instruments.

Use aseptic technique when removing the Cartiva SCI device from the innermost packaging.

Carefully inspect the device and its packaging for any signs of damage, including damage to the sterile barrier. Do not use Cartiva SCI devices if the packaging is damaged or the implant shows signs of damage.

Use care when handling the Cartiva device to ensure that it does not come in contact with objects that could damage the implant. Damaged implants are no longer functionally reliable.

The Cartiva SCI device should not be used with components or instruments from other manufacturers.

Cartiva SCI device should not be re-used or re-implanted. Ensure proper alignment and placement of device components as misalignment may cause excessive wear and/or early failure of the device.